

A new "Molecular Scanner" design for interfacing gel electrophoresis with MALDI-TOF


ThP30 672

Stephen J. Hattan; Kenneth C. Parker; Marvin L. Vestal SimulTof Corporation, Sudbury, MA

Introduction

The concept of a "Molecular Scanner" introduced by Hochstrasser et. al.^{1,2} was intended to provide an efficient mechanism for directly coupled protein gel electrophoresis with detection by MALDI mass spectrometry. Shown here is a new design for the "molecular scanner" incorporating novel membranes that demonstrate efficient capture and digestion of proteins blotted directly for gels. The resulting peptides are eluted from the membranes, detected by MALDI MS and interpreted by peptide mass fingerprinting³ (PMF) with MS/MS confirmation. This work demonstrates the feasibility of this simple analytical workflow for doing top-down proteomics.


Molecular Scanner Concept



- Old**
- PAGE separated protein are orthogonally blotted from gel
 - Proteins pass through an immobilized enzyme membrane
 - Resulting peptides are captured on a hydrophobic membrane
 - Membrane is coated with Matrix and analyzed by MALDI MS
- New**
- new membrane design using silanized glass filter paper
 - efficient protein digestion
 - high capacity peptide capture
 - captured peptide eluted from membrane prior to analysis



Membrane Construction

- Membranes are constructed by the silanization of commercial glass filter paper (Millipore, Pall).
- variety of silanes (Gelest Inc.) were used to create both immobilized enzyme and hydrophobic capture membranes





Membrane Characterization

- Enzyme membrane digestion time
- membrane sampled randomly
 - 3 time points 15 min, 2 hr, 4 hr (BSA sample)
 - 4 samples each time point measured in duplicate



Protein digestion is efficient and reproducible but takes time

- Capture membrane (C8) binding capacity
- membrane sampled randomly
 - 4 sample loads 1, 5, 10, 20 pmol digested BSA onto ~5mm² area
 - 10 sample each load measured in duplicate
 - breakthrough collected and analyzed



5mm² area of membrane consistently binds 20pmol of digested protein with full peptide recovery



Blotting Experiment

- experiment run under normal Western blot conditions
- with and without our capture and digestion membrane

1) yeast lysate loaded in all lanes
 2) clear gel = blotting efficiency
 3) with digestion/capture membranes, protein does not reach PVDF
 4) capture / digestion membrane removed washed and sliced
 5) individual fractions eluted separately (75% Acetonitrile)
 6) all fractions analyzed MALDI-TOF MS (next frame)

Results

34 fractions eluted and analyzed by MALDI-TOF MS



Different fractions contain peptides originating from the proteins in that band --good for PMF--

| Rank | Symbol | Organism | acc. # | AA length | Protein name | # peptides | TrSM | %CM | %IM | ppw |
|------|---------|----------|---------|-----------|-----------------------------|------------|--------|------|-----|-----|
| 0 | TRY_cow | cow | TRY1_BO | 246 | trypsin bovine | 4 | 205137 | 71.7 | 5.5 | 2.4 |
| 1 | RPL35A | Saccharo | P39741 | 120 | 60S ribosomal protein L35 | 5 | 107385 | 31.9 | 8.5 | 0.5 |
| 2 | RPS24A | Saccharo | P26782 | 135 | 40S ribosomal protein S24 | 4 | 57720 | 21.2 | 5.8 | 0.5 |
| 3 | RPL26A | Saccharo | P05743 | 127 | 60S ribosomal protein L26-A | 3 | 36095 | 24.4 | 3.1 | 0.7 |
| 4 | RPS17A | Saccharo | P02407 | 136 | 40S ribosomal protein S17-A | 2 | 31083 | 24.0 | 2.8 | 1.8 |
| 5 | RPL14B | Saccharo | P38754 | 138 | 60S ribosomal protein L14-B | 2 | 21337 | 22.5 | 9.2 | 1.2 |
| 6 | RPS16A | Saccharo | P40213 | 143 | 40S ribosomal protein S16 | 3 | 18959 | 38.9 | 3.4 | 7.7 |
| 7 | RPL4B | Saccharo | P49626 | 362 | 60S ribosomal protein L4-B | 4 | 11957 | 22.4 | 2.6 | 2.2 |
| 8 | RPS19B | Saccharo | P07281 | 144 | 40S ribosomal protein S19-B | 2 | 11695 | 18.7 | 8.6 | 0.8 |
| 9 | RPS14A | Saccharo | P06367 | 137 | 40S ribosomal protein S14-A | 3 | 11211 | 21.0 | 9.3 | 6.8 |

PMF results from fraction 9 (middle above) leads to multiple protein identifications ID's confirmed by MS/MS

Conclusions & Future

Continue to improve membrane design
 Automate membrane fractionation and elution

References

- 1) Binz P-A, et al. (2004) Current Opin in Biotech.; **15**: 17-23
- 2) Nadler TK, et al. (2004) Anal. Biochem.; **332**: 337-348.
- 3) Parker KC (2002) Scoring Methods in MALDI Peptide Mass Fingerprinting.; JASMS: **13**: 22-39.