Comparison of chicken light and dark meat using LC MALDI-TOF mass spectrometry as a model system for biomarker discovery

Jie Du; Stephen J. Hattan ; Kenneth C. Parker

VIC Instruments Corporation, 60 Union Ave, Sudbury, MA

Introduction

Biomarker discovery is challenged by the need to distinguish subtle differences in components contained within complex samples. To date, most proteomic studies of muscle have been performed on human or rodent samples. However, chicken contains muscles distinguishable by eye, *i.d.* light and dark meat, which are expected to contain differences at the protein and protein isoform level Presented here is a comparative study of light and dark meat as a model system for muscle proteomics. Our LC-MALDI analytical platform with high resolution TOF instrumentation (R>40,000) coupled with mass accuracy of <2 ppm across the peptides mass range allows for potential biomarkers to be determined at the MS level of analysis. The protein IDs are confirmed by MS/MS.

Methods

Work Flow

Chicken light and dark meat samples were obtained from a local supermarket. The sample preparation and characterization work flow is illustrated in Fig. 1.



Fig. 1 Work flow chart



Fig. 2 High performance mass spectrometer

Mass Spectrometry

All spectra were acquired on the high performance MALDI-TOF LC MALDI MS: MALDI MS analysis was performed on peptide mixture separated on spectrometers built in our company. ^[1] (Fig. 2). The mass bin was LC with long gradient. **Fig. 4** shows the TIC of 3-hour separation of chicken dark set at 1 ns; average of 1000 laser shots; laser frequency 1 kHz; meat peptides. Fig. 5 gives an example of the MS spectrum showing high resolution scan speed 1 mm/s; mass range 100-2700 Da; operation gas of the instrument. PMF results are summarized in Table 1. pressure 2×10⁻⁸ Torr.

PMF Search Conditions

A peak list consisting of masses and intensities was submitted to our own internal PMF program, ^[2] which iteratively identifies proteins starting from the most abundant. It calculates theoretical intensities from tryptic peptides and decreases protein scores for peptides that are not found. The database consists of ~14,000 sequences and was compiled from UniProt Swiss-Prot and TrEMBL for birds.

Results

Protein Characterization

Protein concentration: UV

absorbance at 280 nm was calibrated using BSA standards. Based on the calibration function, the concentrations of chicken light and dark meat were estimated to be 1.13 mg/mL and 0.62 mg/mL respectively.

Protein digestion efficiency: SDS-PAGE separation of the whole proteins protein samples. Lanes are

before and after trypsin digestion suggests the digestion is efficient and consistent as shown in Fig. 3.



Fig. 3 SDS-PAGE gel of digested and whole chicken numbered from left to right. Lane 1-3 dark meat after digestion, 4-5 dark meat, 6-8 white meat after digestion. 9-10 white meat

Peptide Characterization



Fig. 4 3-hour LC chromatogram for chicken dark meat peptides



Table 1 ID of myosin isoforms by PMF and MS/MS in 3 LC runs each. Number of peptides matched is listed with the protein rank in parentheses. Proteins specific to light or dark meat are color coded. Myosin heavy chain isoforms (MYH) are in pale blue.

R: 48032 Fig. 5 High resolution MS spectrum

R: 49578

R: 45036

of 1476 fraction (R>40,000)

• One dominant MYH in light meat, several alternative MYH in dark meat: PMF occasionally selects different MYHs as the best fit

- 2 MYL isoforms detected
- Green proteins associate with actin. Yellow proteins are general metabolic enzymes in light meat

 Heterogeneity of MYH in dark mea may suppress some proteins.

MS/MS: LC fractions having different peptide masses in MS scan were automatically picked for MS/MS analysis using the software developed by our company. Fig. 6 shows an example of one MS/MS spectrum. Peptide masses that differ in light and dark meat were sequenced and tabulated in Table 2, which confirms the results from PMF.



Conclusion

1382 1

1340

- MS followed by PMF can guickly identify potential biomarkers
- MALDI-TOF MS/MS instrument helps confirm PMF ID
- LC separation increases MS/MS IDs but has little effect on PMF IDs
- Chicken light and dark meat differ by MYH isoforms

Reference

[1] Vestal, M. L. and Hayden, K. Int. J. Mass Spectrom. 2007, 268, 83-92

[2] Parker, K. C.; Garrels, J. I.; Hines, W.; Butler, E. M.; McKee, A. H. Z.; Patterson, D. Martin, S. *Electrophoresis*, 1998, 19(11), 1920-1932

Funding: SBIR Grant 5R44GM079833, 5R44RR025705, 1R43RR025257 5R44GM079832

WP 651

-AA	sequence	isoform	dark	white
759	AGLIGVLEEMR	MYH7	26	
759	AGLLGLEEMR	MYH		56
887	NDLQLQVQAEADSLADAEER	MYH	22	66
887	NDLQLQVQAEADALADAEER	MYH	13	
117	TEELEEEIEAER	MYH6	60	
117	IEELEEEIEAER	MYH	84	75
282	LQTETGEYSR	MYH		72
282	LQTESGEYSR	MYH4	46	
308	qAFTQQIEELK	MYH3	29	
308	qGFTQQIEELK	MYH		34
406	NLQQEIADLTEQIAEGGK	MYH		23
406	NLQQEISDLTEQIAEGGK	MYH	38	
423	LQTEIEDLSVDLER	MYH8	36	
423	LQNEVEDLMVDVER	MYH		61
423	LQNEVEDLMIDVER	MYH3	73	55
680	ANLLQAEIEELR	MYH3	70	
680	ANLLQAETEELR	MYH6	50	
680	ANLLQAEVEELR	MYH	11	63
703	VAEQELMDASER	MYH3	68	
703	VAEQELLDATER	MYH		67
703	VAEQELLDASER	MYH5	25	
703	LAEQELLEATER	MYH8	29	
732	LETDIVQIQSEMEDTIQEAR	MYH		39
732	LETDIAQIQSEMEDTIQEAR	MYH3	93	
732	LESDISQIQSEMEDTIQEAR	MYH5	78	27
732	LETDIVQIQSEMEDTIQEAR	MYH		74
851	ELTYQZEEDR	MYH		70
851	ELTYQSEEDR	MYH4	51	
901	IQHELEEAEER	MYH	60	57
901	VQHELDDAEER	MYH8	20	